

Muon Lifetime Experiment Checkout Procedure

01/19/2005

A guide to a quick and interesting
way to setup the experiment

The Hypothesis

- Random noise in a coincidence experiment with 2 counters can be quantified

$$Rate_{Coincidences}^{Random} = Rate_1 Rate_2 (Width_1 + Width_2)$$


- The sum of the widths of the logic pulses represents the total amount of time over which a coincidence works. (sometimes called the resolving time)
- Imagine the rate of the first counter times this sum of widths as the probability that the coincidence window is open.
- The rate of the second counter times this probability gives you the rate of the accidental coincidences.
- Example: The probability of getting tails on a coin is about $\frac{1}{2}$. If we flip the coin once every 3 seconds, the rate (on average!) we can expect to see tails is: $(1/3s)(1/2) =$ a tail every 6 seconds.

This should work for a single counter too!

- In the muon experiment, we are taking a signal, delaying it, stretching it and putting it in coincidence with itself

Pulse-o-matic diagram of muon set-up

Original Signal from Phototube is turned into a digital signal in the discriminator



Digital pulse is split 3 ways:

To scope




To gate




To
Coincidence



Gate gets delayed and stretched



So if we get a second pulse (indicated by dotted line) the result from the coincidence unit looks like:




The pulse labeled "1" acts like the first counter, the pulse labeled "2" acts like the second counter, and the total width is essentially the stretched width.

Prove it!

- 3 parameters can be varied:
 - 1) The High Voltage on the tube (do not exceed 1600V)
 - 2) The threshold on the discriminator
 - 3) The width of the stretch (Don't change the discriminator width!)
- If the muon experiment is in the “standard set up”
 - 1) The High voltage is 1300V
 - 2) The Discriminator threshold is set to 45mV (reads as 0.450 V on a DVM for easier DVM-ing)
 - 3) The width of the stretch is ~20 microseconds (note: this width must be returned to its original value before data taking starts.
 - 4) There are about 800 single counts/min and 2-10 coincidences/minute (Average is about 6/minute)

What it looks like



Set High Voltage Here

Counts single pulses

Counts Coincidences

Discriminator (Set Threshold here)

Stretcher (Set total Width here)

Setting Threshold with a tiny screwdriver

Measuring Threshold with a DVM



Use the knob on the right (reads 300 now) for HV



Suggested Procedure

- Record the initial setup conditions
- Calculate the expected number of accidental coincidences you expect in a minute with this setup

$$Rate_{Coincidences}^{Random} = Rate_{Single}^2 (Width_{Stretched})$$

- For the “Standard” setup, this is about (1/5) minute
- Turn the knob on the stretcher a couple clicks clockwise or until you notice a substantial increase in the number of coincidences (2 clicks is about a factor of 100)
- Record data for different conditions, varying the 3 parameters from the previous page. Take enough data with enough conditions to verify or refute the above equation
- You’ll notice that coincidences may remain very stable over a large range of conditions. That’s probably signal. I.e.

$$Rate_{Coincidences} = Rate_{Coincidences}^{Random} + Rate_{Coincidences}^{Signal}$$

- Is the “standard” setup still a good choice?